



CFC FREE Cold Chain Equipments

Truly World Class



LONG RANGE

**COLD BOX
BRAND : APEXIN
WHO PQS NO.: E004/010
MODEL : AICB-444L**



SHORT RANGE

**COLD BOX
BRAND : APEXIN
WHO PQS.: E004/030
MODEL : AICB-243S**

SPECIFICATIONS

MODEL AICB-444L

MODEL AICB-243S

	Long	Short
Range	Long	Short
Storage Volume	18 L	8 L
Weight Fully Loaded	48.87 kg	22.41 kg
Weight Empty	19.93 kg	11.57 kg
External Dimensions (LxWxH)	76.1 x 61.1 x 51.3 cm	54.4 x 44.5 x 42.0 cm
Internal Dimensions (LxWxH)	52 x 37 x 28.2	34.6 x 28.0 x 24.0 cm
Storage Dimensions (LxWxH)	44.8 x 30.0 x 16.7 cm	27.0 x 20.4 x 15.0 cm
Lid Type & Fixing	Removable lid	Fixed Hinges
External material	Polyethylene	Polyethylene
Internal Lining Material	Polyethylene	Polyethylene
Insulation material	Polyurethane	Polyurethane
Insulation Thickness	100 mm	65 mm
Number of Ice Packs	44 x 0.4 L	24 x 0.3 L
Cold life without opening at 43°C	140 hours	84 hours 7 mn

Distributed by:



PT. SIGMA BIMED

Indonesia Sole Distributor

Jl. Agung Utara Raya Blok A 36D/42 Sunter Agung Jakarta Utara 14350 Indonesia

Telp. (62-21) 64716568, 87960123, 87960456

Fax. (62-21) 64716566, 87961123, 87961456

Email : sigmabimed@indo.net.id, sales@sigmabimed.com

Website : www.sigmabimed.com