



CFC FREE Cold Chain Equipments

Truly World Class

**VACCINE CARRIER
MODEL : AIDVC-24
SHORT RANGE
BRAND : APEXIN
WHO PQS NO.: E004/011**



**VACCINE CARRIER
MODEL : AIVC-44
LONG RANGE
BRAND : APEXIN
WHO PQS NO.: E004/012**



SPECIFICATIONS

	MODEL AIDVC-24	MODEL AIVC-44
Storage Volume	0.9 L	1.6 L
Weight Fully Loaded	2.27 kg	4.33 kg
Weight Empty	1.31 kg	1.93 kg
External Dimensions (LxWxH)	25.0 x 18.0 x 21.4 cm	25.0 x 25.0 x 30.0 cm
Internal Dimensions (LxWxH)	16.5 x 10.0 x 12.0 cm	9.96 x 9.96 x 16.7 cm
Storage Dimensions (LxWxH)	16.5 x 10.0 x 5.0 cm	9.0 x 9.0 x 16.7 cm
Lid Type & Fixing	Removable lid	Removable lid
External material	Polyethylene	HDPE (High Density Polyethylene)
Internal Lining Material	Polyethylene	HIPS (High Impact Polystyrene)
Insulation material	Polyurethane	Polyurethane
Insulation Thickness	32 mm	35 mm
Number of Ice Packs	2 x 0.4 L	4 x 0.4 L
Cold life without opening, tested by WHO approved laboratory	21 hours 30 mn at 43°C	39 hours at 43°C

Distributed by:



PT. SIGMA BIMED

Indonesia Sole Distributor

Jl. Agung Utara Raya Blok A 36D/42 Sunter Agung Jakarta Utara 14350 Indonesia

Telp. (62-21) 64716568, 87960123, 87960456

Fax. (62-21) 64716566, 87961123, 87961456

Email : sigmabimed@indo.net.id, sales@sigmabimed.com

Website : www.sigmabimed.com