



Japanese Medical Device Manufacturers
Certification Number : 27B2X00214

COOL ALASTAR®

Cold Bench



Fridge type

Temperature range:
2 ~ 8°C

Model : DCA-B120R
ISO : 9001
ISO : 13485

Features

▶ Giving enough time to process . . .

- Safety control device in the Fridge mode gives alarm at 2°C ,shuts down the compressor to prevent temperatures lower than about 0°C on the table as super cooling
- One-touch operation is easy to use
- Electronic digital thermostat and temperature indicator enables reading easy
- Special cooling sensor controls the frost automatically to maintain a constant cooling performance
- Attached casters enable to smooth movement
- 2 independent digital alarms are designed to detect abnormality


Use

▶ The low temperature table has been developed for the routine work of . .


- Quality control of Blood products
- Temporary storage of Vaccine
- Clinical research
- Biological research

Dimensions

DCA - B60R



DCA - B120R



Specification

Model	Fridge type			Fridge-Freeze switch type	
	DCA-B60R	DCA-B120R	DCA-B150R	DCA-B120RF	DCA-B150RF
certificate No.	21900BZX00014000				
Temperature range	2~8°C			below- 20°C or 2~8°C	
Outside dimensions (mm)	W 600	W 1200	W 1500	W 1200	W 1500
	D 750	D 750	D 750	D 750	D 750
	H 850-860	H 850-860	H 850-860	H 850-860	H 850-860
Materials	TABLE: Stainless steel SUS-304 Exterior: Colored Bonderizing zinc steel				
Refrigerant Gas	HFC134a (CFC Free)		HFC404A (CFC Free)		
Cooling method	Forced circulation blow-up				
Control method	Electronic digital thermostat, microcontroller				
Alarm device	2 digital Thermo-alarm devices, High temp +8°C ; Low temp +2°C				
Safety device	Shut down compressor when overcooling			Shut down compressor when overcooling in fridge mode	
Electricity	1 phase 220V 50Hz with transformer				

Attention

- Please shut off the power supply after working to defrost
- Cleaning the condenser filter once a month
- Inspection the equipment and alarm temperature once a year
- Recommendation RH: 40-70% , AT: 10-28°C
- Read the instruction manual before using
- Made in Japan

◆ Manufacturer

 **Daido Industries INC.**

◆ Distributor

 **PT. SIGMA BIMED**

Indonesia Sole Distributor
 Jl. Agung Utara Raya Blok A 36D/42
 Sunter Agung, Jakarta 14350, Indonesia
 Telp : (62-21) 64716568, 70240088, 70240099
 Fax : (62-21) 64716566, 87961123, 87961456
 Email : sigmabimed@indo.net.id | sigmaeng@indo.net.id
 Website : http://www.sigmabimed.com
 Kemenkes RI AKL : 20209411199