

Preservation above-zero

Blood Bank Refrigerators

Specially designed to provide precise and reliable temperature control that meets various criteria for whole blood and blood derivatives.



MBR-506D(H)

\$ " ! % # \$ " ! & ' () ' * *
 ! +
 % # \$ %
 &
 ! #
 , ! " # " - / \$ %
 *
 0


Alarm and Safety Functions

To ensure the safety of precious blood supply, the Panasonic MBR Series provides the following safety functions:

- **Audible and flashing LED visual alarms** with remote alarm contacts, in case of power failure, high or low temperature condition, or any thermal sensor abnormality.
- Door alarm and key lock are standard features.
- Re-activating buzzer, lamp and remote alarm contact. (30 min. after buzzer stops)
- Built-in temperature recorder.

Environmentally Friendly

Panasonic refrigerators feature commercially-available HCFC-free, HFC refrigerants and HCFC-free insulation.



MBR-506D(H)

Model No.	MBR-506D(H)
Exterior dimensions (W x D x H)	800 x 832 x 1810 (mm) / 31.5 x 32.8x71.3
Net weight	185kg
Storage capacity	120 bags (450ml), with MBR-56B
Interior cabinet	Stainless steel
Outer door	Triple layer glass windows
Inner door (acrylic)	5
Casters	4
Refrigerant	HFC
Air circulation	Forced air circulation
Temperature control	Microprocessor
Power failure alarm	Audible and visual alarm (9 hours), Automatic rechargeable battery
Remote alarm contact	DC24V, 1A Temperature alarm or power failure alarm
Access port	1
Options	Recording paper (RP-06-PW) Basket (MBR-56B-PW) Max. 20pcs/unit 450ml bags x 6pcs/basket

Panasonic Healthcare Co., Ltd., Gunma Factory is certified for:
 Quality management system: ISO9001
 Medical devices quality management system: ISO13485

Panasonic Healthcare Co., Ltd., Gunma Factory is certified for:
 Environmental management system: ISO14001

DISTRIBUTED BY: