

Panasonic

CHEST FREEZER SF-PC1497 (Hard Top Lid)



FEATURES

- Solid door chest freezer
- Inside temperature -20°C at load line
- The insulation is made from polyurethane CFC-Free, thickness 75 mm
- Temperature control is equipped at the front for convenient using
- Coil is made from copper is durable and anticorrosion, long life using

Temperature range	: -20°C
Capacity	: 382 liter
Temperature control	: Thermostat
Electrical	: 220V/50Hz
Compressor (Watt)	: 200
Refrigerant	: R-134 (CFC-Free)
Internal dimensions (WxDxH)	: 1180x500x698 mm
External dimensions (WxDxH)	: 1335x747x945 mm
Net weight	: 74 Kg

IMPORTED AND DISTRIBUTED BY :



PT. SIGMA BIMED

Indonesia Sole Distributor

Jl. Agung Utara Raya Blok A 36/42 Sunter Agung Jakarta Utara 14350 Indonesia

Telp. (62-21) 64716568, 70240088, 70240099, 87960123, 87960456

Fax. (62-21) 64716566, 87961123, 87961456

Email : sigmabimed@indo.net.id, sales@sigmabimed.com

Website : www.sigmabimed.com